

CytoSure



A Sysmex Group Company

Application Note

The use of the InnoScan 710 scanner and Mapix software with CytoSure microarrays

Adrianna Lagraulet (Innopsys)

The use of the InnoScan[®] 710 scanner and Mapix[®] software with CytoSure microarrays

Adrianna Lagraulet (Innopsys)

Introduction

OGT's CytoSure[®] products provide researchers with the complete solution for highly accurate detection and interpretation of chromosomal abnormalities. They facilitate high-throughput analysis of the whole human genome to accurately identify the genetic variation underlying hundreds of genetic disorders and congenital syndromes that cause developmental delay. Critical to success, is the use a scanner which offers ease-of-use, high precision and reliability of results with a wide range of array formats and densities.

The InnoScan 710 scanner from Innopys is a widely used, low cost, portable microarray scanner capable of scanning microarray slides with a resolution of 3 μm . It is the fastest scanner on the market with an adjustable reading speed from 10 to 35 lines per second. Offering reliable and traceable results, with a barcode reader, InnoScan 710 provides accurate control of acquired images and is supplied with a validation slide and associated Mapix[®] software that verifies scanner performance as well as offering a sophisticated feature extraction function.

With two high performance lasers with an excitation of 532nm and 635nm, the InnoScan 710 scanner complements OGT's CytoSure range of microarray slides, reagents and analysis software. This application note illustrates the use of the InnoScan 710 scanner and Mapix with a range of OGT CytoSure cytogenetic array formats.

Materials and Methods

A range of OGT CytoSure slides, the CytoSure ISCA v2 (8x60k and 4x180k) arrays, the CytoSure Embryo Screen array (8x60k) and the CytoSure Medical Research Exome array (1x1M), were processed according to OGT protocols (ogt.com). Briefly, in the case of the CytoSure Embryo Screen Arrays, DNA was first amplified using the PicoPLEX™ kit (Rubicon). All other slides used 1µg of starting DNA. The DNA was then labelled using the CytoSure Genomic DNA Labelling Kit and purified with the supplied CytoSure columns. Sample DNA was labelled with Cy3 and Reference DNA with Cy5. After purification, the DNAs were pooled and dried in a SpeedVac™ until dry. The DNA was resuspended in water, CotI, 5x blocking agent (Agilent) and 2x High RPM hybridization buffer (Agilent) and hybridised to OGT's microarray slides.

After hybridisation, the slides were washed and scanned with the InnoScan 710 microarray scanner. Innopsys has developed an auto-setting function dedicated to oligo slides, either 8x60k or 4x180k. Using these slide configurations, all of the scan parameters are automatically adapted to slide format to get optimal signals in compliance with OGT's QC metrics. The Oligo 8x60k and Oligo 4x180k scan configurations are available from Mapix version 7.4.0.

Following scanning, the image was feature extracted using the Mapix software. This was carried out by importing the GAL file supplied with the CytoSure array and feature extracting all the arrays on the slide simultaneously. The resulting .gpr files generated were then imported into CytoSure Interpret Software. Normalisation was carried out in the software using the Segmentation-based normalisation. The process is summarised in Figure 1.




Figure 1: Summary of the process used in this study.

Results

CytoSure ISCA v2 8x60k array

The CytoSure ISCA v2 8x60k array comprises a research-validated collection of specific probes that enable reliable detection of copy number (CN) changes for a variety of genetic disorders. They offer excellent coverage with probes spaced throughout the genome and complete coverage, with increased probe density, in areas recommended by ClinGen (formerly the ISCA consortium). Thus ensuring uncompromising detection of CN changes in these cytogenetically important regions.

A CytoSure ISCA v2 8x60k slide was hybridised with DNA samples with a wide range of sized chromosomal gains and losses, ranging from large, such as NA00072 with a ~25Mb deletion, to small intragenic deletions. The results are shown in the Table below and in Figure 2.

Sample	Aberration	DLRs	Green Signal	Red Signal	Detection
NA22624	arr 11p12p11.2(40433344-46031324)x1	0.15	3099	2292	Yes
NA18324	arr 17p11.2(16697859-20294610)x1	0.16	2472	3300	Yes
NA13434	Exon 3-4 PLP1	0.17	2576	3194	Yes
NA14117	arr 5p15.33p14.3(68519-22367289)x1	0.17	2372	3093	Yes
NA06937	10pter>10p1	0.17	3057	2833	Yes
NA00072	arr 4p16.3p15.2(1618261-25512740)x1	0.16	2744	3253	Yes

Table 1: Results from the CytoSure ISCA v2 8x60k array.




Figure 2: Detection of a 3.6Mb deletion on chromosome 17 in sample NA18324.

CytoSure ISCA v2 4x180k array

The CytoSure 4x180k ISCA v2 array is similar to the 60k array in that it offers excellent coverage with probes spaced throughout the genome and complete coverage, with increased probe density, in areas recommended by ClinGen. The larger probe set of ~180k probes, offers superior CN resolution of 1 probe every 19kb in targeted regions and 1 probe every 25kb along the backbone, (compared to 1 probe every 40kb and 70kb, respectively, with the 8x60k array). The performance of the InnoScan710 when scanning CytoSure ISCA v2 4x180k slides was assessed. All arrays gave good or excellent DLRs (Table 2).

Sample	DLRs	Green Signal	Red Signal
Male / Female	0.21	1727	1371
Male / Female	0.18	1555	1638
Male / Female	0.20	2166	1924
Male / Female	0.23	1698	1555

Table 2: CytoSure ISCA v2 4x180k results, arrays have been scanned with an InnoScan 710 scanner and feature extracted using Mapix.

CytoSure Embryo Screen 8x60k array

The CytoSure Embryo Screen array utilises long oligo array comparative genomic hybridisation (aCGH) for superior aneuploidy and copy number detection. The array content has been optimised to work with small amounts of DNA amplified from a single cell, specifically single cells isolated from a pre-implantation embryo, following amplification by the Rubicon PicoPLEX kit. In this experiment 50pg of DNA was amplified using the PicoPLEX kit. The amplified DNA was then labelled and hybridised as detailed in the CytoSure Embryo Screen protocol. The slide was scanned using the InnoScan 710 scanner, feature extracted using Mapix and imported into CytoSure Interpret Software. Analysis was carried out using the PGS interpretation module and the results shown in Table 3 and Figure 3. All arrays gave good Signal to Noise Ratios (SNRs), while the DLRs were well below 1.0 which is classed as ‘excellent’ for PicoPLEX amplified material and is sufficient to detect aberrations of 15Mb and above.

Sample	Aberration	DLRs	Green Signal	Red Signal	Detection
Male / Female reference	n/a	0.60	4312	5487	Yes
NA14117	arr 5p15.33p14.3(68519-22367289)x1	0.63	4555	3828	Yes
NA00072	arr 4p16.3p15.2(1618261-25512740)x1	0.64	3716	5740	Yes (manual)
NA01921	arr 21q11.2q22.3(13583116-46921373)x3	0.73	4508	4463	Yes

Table 3: CytoSure Embryo Screen Array results




Figure 3: Detection of a 15Mb deletion on chromosome 5 in sample NA14117, using a CytoSure Embryo Screen Array.

CytoSure Medical Research Exome array (1x1M)

The CytoSure Medical Research Exome array is a highly targeted exon-focussed array capable of detecting medically relevant microdeletions and microduplications. The array makes an ideal complement to an exome sequencing approach to provide a comprehensive mutation spectrum analysis in rare disease.

NA13434 DNA was labelled and hybridised to the array using the standard protocol and scanned using the InnoScan 710 scanner. Results are shown in Table 4 and Figure 4.

Sample	Aberration	DLRs	Green Signal	Red Signal	Detection
NA13434	PLP1 deletion in exon 3 and 4	0.225	1785	1006	Yes

Table 4: CytoSure Medical Research Exome array result




Figure 4: Detection of an aberration within the PLP1 gene using the CytoSure Medical Research Exome array scanned with the InnoScan 710 scanner.

Conclusions

This application note illustrates how the InnoScan 710 scanner can be used effectively to scan a wide range of OGT's CytoSure arrays to accurately detect chromosomal copy number variation. The Mapix software feature extracts the scanned image and can be coupled with the CytoSure Interpret Software to analyse the results in a fast, efficient process. It is recommended to use the oligo 8x60k or oligo 4x180k dedicated scanner auto settings, adapted to CytoSure 8x60k and 4x180k slides, respectively.

Ordering information

UK +44 (0) 1865 856800

US +1 914 467 5285

contact@ogt.com

ogt.com



A Sysmex Group Company

**What binds us,
makes us.**

Oxford Gene Technology Ltd., Begbroke Science Park, Woodstock Road, Begbroke, Oxfordshire, OX5 1PF, UK

CytoSure: This product is provided under an agreement between Agilent Technologies, Inc. and OGT. The manufacture, use, sale or import of this product may be subject to one or more of U.S. patents, pending applications, and corresponding international equivalents, owned by Agilent Technologies, Inc. The purchaser has the non-transferable right to use and consume the product for RESEARCH USE ONLY AND NOT for DIAGNOSTICS PROCEDURES. It is not intended for use, and should not be used for the diagnosis, prevention, monitoring, treatment or alleviation of any disease or condition, or for the investigation of any physiological process, in any identifiable human, or for any other medical purpose. OGT™ and CytoSure® are trademarks of Oxford Gene Technology IP Ltd. InnoScan® and Mapix® are registered trademarks of Innopsys. This document and its contents are © Oxford Gene Technology IP Limited – 2021. All rights reserved.