



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of: Hazardous Substances (Safety Data Sheets) Notice 2017. This notice is issued by the Environmental Protection Authority under sections 75 and 76(1)(b), (f), (g) and (h) of the Hazardous Substances and New Organisms Act 1996

Issuing Date 15-Oct-2021

Revision date 04-Dec-2024

Revision Number 2

## Section 1: Identification

### Product identifier

Product Name Rubber Solution Glue

Product Code(s) PCA005

### Other means of identification

### Recommended use of the chemical and restrictions on use

Recommended use Adhesives

Uses advised against None known

### Details of the supplier of the safety data sheet

#### Importer

Sysmex New Zealand Limited  
Level 3, 103 Carlton Gore Rd  
New Market  
Auckland 1023, New Zealand  
+64-9-630-3554/ 0800797639

#### Manufacturer

CytoCell Ltd., Oxford Gene Technology  
418 Cambridge Science Park, Milton Road,  
Cambridge  
CB4 0PZ, United Kingdom  
T: +44 (0)1223 294048  
F: +44 (0)1223 294986  
probes@cytoCell.com  
<http://www.ogt.com>

E-mail address regulatory@sysmex.co.nz

### Emergency telephone number

Emergency telephone For Sysmex Supply Chain support or Product Related Enquiries: +64 9 6303554 / 0800797639 (Mon to Fri – 8.30 am to 5.00 pm)  
For any spillage or clean up issues: CHEMCALL 0800 243 622 (24 hours – 365 days)  
National Poison Centre 0800 764 766 (0800 POISON)

## Section 2: Hazard identification

### Classification of the substance or mixture

Flammable liquids	Category 2
Aspiration hazard	Category 1
Skin corrosion/irritation	Category 2
Specific target organ toxicity (single exposure)	Category 3
Category 3 Target organ effects: Narcotic effects.	
Hazardous to the aquatic environment - acute	Category 1
Hazardous to the aquatic environment - chronic	Category 1

### Label elements

**Signal word**

DANGER

**Hazard statements**

Highly flammable liquid and vapour  
 Causes skin irritation  
 May cause drowsiness or dizziness  
 Very toxic to aquatic life with long lasting effects  
 May be fatal if swallowed and enters airways

**Precautionary Statements - Prevention**

Use explosion-proof electrical/ ventilating/ lighting/ equipment  
 Wash face, hands and any exposed skin thoroughly after handling  
 Avoid breathing dust/fume/gas/mist/vapours/spray  
 Use only outdoors or in a well-ventilated area  
 Avoid release to the environment  
 Ground and bond container and receiving equipment  
 Use non-sparking tools  
 Take action to prevent static discharges  
 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking  
 Keep container tightly closed  
 Wear protective gloves/clothing and eye/face protection  
 Keep cool

**Precautionary Statements - Response**

Specific treatment (see supplemental first aid instructions on this label)  
 If skin irritation occurs: Get medical advice/attention  
 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]  
 Wash contaminated clothing before re-use  
 IF INHALED: Remove person to fresh air and keep comfortable for breathing  
 Call a POISONS INFORMATION CENTRE or doctor if you feel unwell  
 IF SWALLOWED: Immediately call a POISONS INFORMATION CENTRE or doctor  
 Do NOT induce vomiting  
 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish  
 Collect spillage

**Precautionary Statements - Storage**

Store locked up  
 Store in a well-ventilated place. Keep container tightly closed

**Precautionary Statements - Disposal**

Dispose of contents/container to an approved waste disposal plant

**Other hazards which do not result in classification**

The product does not contain any substance(s) classified as PBT or vPvB. May be harmful in contact with skin.

**Section 3: Composition/information on ingredients**

Chemical name	CAS No.	Weight-%
n-Heptane	142-82-5	60 - 90
Ethanol	64-17-5	10 - 40

## Section 4: First-aid measures

### Description of first aid measures

<b>General advice</b>	Show this safety data sheet to the doctor in attendance Immediate medical attention is required
<b>Inhalation</b>	Remove to fresh air. Aspiration into lungs can produce severe lung damage. If breathing has stopped, give artificial respiration. Get medical attention immediately. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. If breathing is difficult, (trained personnel should) give oxygen. Delayed pulmonary oedema may occur.
<b>Eye contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Keep eye wide open while rinsing. Do not rub affected area.
<b>Skin contact</b>	Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. Get medical attention if irritation develops and persists.
<b>Ingestion</b>	Do NOT induce vomiting. Rinse mouth. Never give anything by mouth to an unconscious person. ASPIRATION HAZARD IF SWALLOWED - CAN ENTER LUNGS AND CAUSE DAMAGE. If vomiting occurs spontaneously, keep head below hips to prevent aspiration. Get immediate medical attention.
<b>Self-protection of the first aider</b>	Remove all sources of ignition. Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. Use personal protective equipment as required. See section 8 for more information. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. Avoid contact with skin, eyes or clothing.

### Most important symptoms and effects, both acute and delayed

<b>Symptoms</b>	Difficulty in breathing. Coughing and/or wheezing. Dizziness. Inhalation of high vapour concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting.
<b>Effects of Exposure</b>	No information available.

### Indication of any immediate medical attention and special treatment needed

<b>Note to doctors</b>	Because of the danger of aspiration, emesis or gastric lavage should not be employed unless the risk is justified by the presence of additional toxic substances.
------------------------	---

## Section 5: Fire-fighting measures

**Hazchem code** •3YE

### Suitable Extinguishing Media

**Suitable Extinguishing Media** Dry chemical. Carbon dioxide (CO<sub>2</sub>). Water spray. Alcohol-resistant foam.

**Unsuitable extinguishing media** Do not use a solid water stream as it may scatter and spread fire.

### Specific hazards arising from the chemical

**Specific hazards arising from the chemical** Highly flammable liquid and vapour. Vapours are heavier than air and may spread along floors. Risk of ignition. Keep product and empty container away from heat and sources of ignition. In the event of fire, cool tanks with water spray. Fire residues and contaminated fire extinguishing water must be disposed of in accordance with local regulations.

**Hazardous combustion products** Carbon oxides.

**Special protective actions for firefighters**

**Special protective equipment and precautions for fire-fighters** Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

**Section 6: Accidental release measures**

**Personal precautions, protective equipment and emergency procedures**

**Personal precautions** Evacuate personnel to safe areas. Use personal protective equipment as required. See section 8 for more information. Avoid contact with skin, eyes or clothing. Ensure adequate ventilation. Keep people away from and upwind of spill/leak. ELIMINATE all ignition sources (no smoking, flares, sparks or flames in immediate area). Pay attention to flashback. Take precautionary measures against static discharges. All equipment used when handling the product must be grounded. Do not touch or walk through spilled material.

**Other information** Ventilate the area. Refer to protective measures listed in Sections 7 and 8.

**For emergency responders** Use personal protection recommended in Section 8.

**Environmental precautions**

**Environmental precautions** Refer to protective measures listed in Sections 7 and 8. Prevent further leakage or spillage if safe to do so. Prevent product from entering drains.

**Methods and material for containment and cleaning up**

**Methods for containment** Stop leak if you can do it without risk. Do not touch or walk through spilled material. A vapour suppressing foam may be used to reduce vapours. Dam far ahead of spill to collect run-off water. Keep out of drains, sewers, ditches and waterways. Absorb with earth, sand or other non-combustible material and transfer to containers for later disposal.

**Methods for cleaning up** Take precautionary measures against static discharges. Dam up. Soak up with inert absorbent material. Pick up and transfer to properly labelled containers.

**Precautions to prevent secondary hazards**

**Prevention of secondary hazards** Clean contaminated objects and areas thoroughly observing environmental regulations.

**Section 7: Handling and storage**

**Precautions for safe handling**

**Advice on safe handling** Contaminated work clothing should not be allowed out of the workplace. Regular cleaning of equipment, work area and clothing is recommended. Use personal protection equipment. Avoid breathing vapours or mists. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Use grounding and bonding connection when transferring this material to prevent static discharge, fire or explosion. Use with local exhaust ventilation. Use spark-proof tools and explosion-proof equipment. Keep in an area equipped with sprinklers. Use according to package label instructions. Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Take off contaminated clothing and wash before re-use. In case of insufficient ventilation, wear suitable respiratory equipment.

**General hygiene considerations** Do not eat, drink or smoke when using this product. Contaminated work clothing should not

be allowed out of the workplace. Regular cleaning of equipment, work area and clothing is recommended. Wash hands before breaks and immediately after handling the product. Wear suitable gloves and eye/face protection.

**Conditions for safe storage, including any incompatibilities**

**Storage Conditions** Keep containers tightly closed in a dry, cool and well-ventilated place. Keep away from heat, sparks, flame and other sources of ignition (i.e. pilot lights, electric motors and static electricity). Keep in properly labelled containers. Do not store near combustible materials. Keep in an area equipped with sprinklers. Store in accordance with the particular national regulations. Store in accordance with local regulations. Store locked up. Keep out of the reach of children. Store away from other materials.

**Incompatible materials** Strong acids, Strong bases, Strong oxidising agents.

**Section 8: Exposure controls/personal protection**

**Control parameters**

**Exposure Limits**

Chemical name	New Zealand	Australia	ACGIH TLV	United Kingdom
n-Heptane 142-82-5	TWA: 400 ppm TWA: 1640 mg/m <sup>3</sup> STEL: 500 ppm STEL: 2050 mg/m <sup>3</sup>	TWA: 400 ppm TWA: 1640 mg/m <sup>3</sup> STEL: 500 ppm STEL: 2050 mg/m <sup>3</sup>	TWA: 400 ppm STEL: 500 ppm	TWA: 500 ppm TWA: 2085 mg/m <sup>3</sup> STEL: 1500 ppm STEL: 6255 mg/m <sup>3</sup>
Ethanol 64-17-5	TWA: 200 ppm TWA: 380 mg/m <sup>3</sup> STEL: 800 ppm STEL: 1520 mg/m <sup>3</sup>	TWA: 1000 ppm TWA: 1880 mg/m <sup>3</sup>	STEL: 1000 ppm	TWA: 1000 ppm TWA: 1920 mg/m <sup>3</sup> STEL: 3000 ppm STEL: 5760 mg/m <sup>3</sup>

**Biological occupational exposure limits** This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

**Appropriate engineering controls**

**Engineering controls** Showers  
Eyewash stations  
Ventilation systems.

**Individual protection measures, such as personal protective equipment**

**Eye/face protection** Tight sealing safety goggles.

**Hand protection** Ensure that the breakthrough time of the glove material is not exceeded. Refer to glove supplier for information on breakthrough time for specific gloves. Wear suitable gloves. Impervious gloves.

**Skin and body protection** Wear suitable protective clothing. Long sleeved clothing. Chemical resistant apron. Antistatic boots.

**Respiratory protection** No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

**Environmental exposure controls** Avoid release to the environment.

## Section 9: Physical and chemical properties

### Information on basic physical and chemical properties

#### Appearance

Physical state	Liquid
Colour	Yellowish
Odour	Petroleum
Odour threshold	No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH		No data available
pH (as aqueous solution)		No data available
Melting point / freezing point		No data available
Initial boiling point and boiling range		No data available
Flash point	-4 °C	CC (closed cup)
Flammability		No data available
Flammability Limit in Air		
Upper flammability or explosive limits	6.7 %	
Lower flammability or explosive limits	1.1 %	
Vapour pressure	98760 mmHg	
Relative vapour density		No data available
Relative density	0.722	
Bulk density		No data available
Liquid Density		No data available
Solubility(ies)		No data available
Water solubility	Immiscible in water	
Partition Coefficient (n-octanol/water)		No data available
Auto-ignition temperature		No data available
Decomposition temperature		No data available
SADT (°C)		No data available
Kinematic viscosity		No data available
Dynamic viscosity	400 - 600 cP	@ 20 °C
Particle characteristics		
Particle Size		No data available
Particle Size Distribution		No data available
Explosive properties	No information available.	
Oxidising properties	No information available.	

#### Other information

Molecular weight	No information available
VOC content	660 g/L 660
Softening point	No information available

### Information with regard to physical hazard classes

#### Explosives

Explosive properties No information available.

#### Oxidising properties

No information available.

## Section 10: Stability and reactivity

**Reactivity**

**Reactivity** None under normal use conditions.

**Chemical stability**

**Stability** Stable under normal conditions.

**Explosion data**

**Sensitivity to mechanical impact** None.

**Sensitivity to static discharge** Yes.

**Possibility of hazardous reactions**

**Possibility of hazardous reactions** None under normal processing.

**Conditions to avoid**

**Conditions to avoid** Protect from direct sunlight. Heat, flames and sparks.

**Incompatible materials**

**Incompatible materials** Strong acids, Strong bases, Strong oxidising agents.

**Hazardous decomposition products**

**Hazardous decomposition products** Carbon oxides.

**Section 11: Toxicological information****Acute toxicity****Information on likely routes of exposure****Product Information**

**Inhalation** Specific test data for the substance or mixture is not available. Aspiration into lungs can produce severe lung damage. May cause pulmonary oedema. Pulmonary oedema can be fatal. May cause irritation of respiratory tract. May cause drowsiness or dizziness.

**Eye contact** Specific test data for the substance or mixture is not available. May cause irritation.

**Skin contact** Repeated exposure may cause skin dryness or cracking. Specific test data for the substance or mixture is not available. Causes skin irritation. (based on components).

**Ingestion** Specific test data for the substance or mixture is not available. Potential for aspiration if swallowed. May cause lung damage if swallowed. Aspiration may cause pulmonary oedema and pneumonitis. May be fatal if swallowed and enters airways. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhoea.

**Symptoms** Difficulty in breathing. Coughing and/or wheezing. Dizziness. Redness. May cause redness and tearing of the eyes. Inhalation of high vapour concentrations may cause symptoms like headache, dizziness, tiredness, nausea and vomiting.

**Acute toxicity** No information available.

**Numerical measures of toxicity**

The following values are calculated based on chapter 3.1 of the GHS document:

ATEmix (dermal) 3,333.30 mg/kg

**Component Information**

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
n-Heptane	-	= 3000 mg/kg ( Rabbit )	> 29.29 mg/L ( Rat ) 4 h
Ethanol	= 7060 mg/kg ( Rat )	-	= 116.9 mg/L ( Rat ) 4 h = 133.8 mg/L ( Rat ) 4 h

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

**Skin corrosion/irritation** Classification based on data available for ingredients. Causes skin irritation.

**Serious eye damage/eye irritation** No information available.

**Respiratory or skin sensitisation** No information available.

**Germ cell mutagenicity** No information available.

**Carcinogenicity**

The table below indicates whether each agency has listed any ingredient as a carcinogen.

Chemical name	New Zealand	IARC
Ethanol - 64-17-5		Group 1

**Legend**

**IARC (International Agency for Research on Cancer)**

Group 1 - Carcinogenic to Humans

**Reproductive toxicity** No information available.

**STOT - single exposure** May cause drowsiness or dizziness.

**STOT - repeated exposure** No information available.

**Aspiration hazard** May be fatal if swallowed and enters airways.

**Data used to identify the health effects** Refer to Section 16 for Key literature references and sources for data used to compile the SDS.

**Section 12: Ecological information****Ecotoxicity**

**Ecotoxicity** Very toxic to aquatic life with long lasting effects.

**Aquatic ecotoxicity**

Chemical name	Algae/aquatic plants	Fish	Crustacea
n-Heptane	-	LC50: =375.0mg/L (96h, Cichlid fish)	-
Ethanol	-	LC50: 12.0 - 16.0mL/L (96h, Oncorhynchus mykiss) LC50: >100mg/L (96h, Pimephales promelas)	LC50: 9268 - 14221mg/L (48h, Daphnia magna) EC50: =2mg/L (48h, Daphnia magna)



		LC50: 13400 - 15100mg/L (96h, Pimephales promelas)	
--	--	--	--

**Terrestrial ecotoxicity**

Chemical name	Earthworm	Avian	Honeybees
Ethanol	Acute Toxicity: LC50 0.1 - 1 mg/cm2 (Eisenia foetida, 48 h filter paper)	-	-

**Persistence and degradability**

**Persistence and degradability**            The product is expected to be slowly biodegradable.

**Bioaccumulative potential**

**Bioaccumulation**

**Component Information**

Chemical name	Partition coefficient
n-Heptane	4.66
Ethanol	-0.35

**Mobility in soil**

**Mobility in soil**                            No information available.

**Other adverse effects**

No information available.

**Section 13: Disposal considerations**

**Disposal methods**

**Waste from residues/unused products**

Dispose of product in packaging in a way that is consistent with the EPA Consolidation 30 April 2021 of the Hazardous Substances (Disposal) Notice 2017 and the Act. Treat the substance using a method that changes the characteristics or composition of the substance so that the substance is no longer a hazardous substance; or export the substance from New Zealand as waste. Flammable substances - may not be disposed of into or onto a landfill or sewage facility.

They may only be burnt in certain situations.

Flammable gases, liquids and solids may only be discharged into the environment or landfill as waste if the substance will not at any time come into contact with any explosives, oxidising gases, liquids or solids or organic peroxides; and there will be no ignition source in the vicinity of the disposal site at any time and if the substance were to ignite, no person, or place where a person may legally be, would be exposed to an unsafe level of heat radiation. Substances which are hazardous to human health or corrosive to metals – may be discharged into the environment if a tolerable exposure limit has been set for the substance (or a component of that substance); and the discharge does not, after reasonable mixing, result in the concentration of the substance in an environmental medium exceeding the tolerable exposure limit. If there is no tolerable exposure limit for the substance, then it may only be discharged into the environment if the substance is very rapidly converted to substances that are not hazardous substances. Environmentally

hazardous substances – if the substance, or if it contains a component that is hazardous to the aquatic environment or bioaccumulative and not rapidly degradable, then any component that is bioaccumulative and not rapidly degradable must be removed. The product may only be discharged into the environment if an environmental exposure limit has been set for the substance (or a component of the substance); and the discharge does not, after reasonable mixing, result in the concentration of the substance in an environmental medium exceeding the environmental exposure limit.

**Contaminated packaging**

For packages that have been in direct contact with hazardous substances, the person must ensure that the package is rendered incapable of containing any substance. It must be disposed of in a manner that is consistent with the requirements for disposal of the substance that it contained, taking into account the material the package is manufactured from. Packages may only be reused or recycled if:

- the substance has a physical hazard other than corrosive to metal, and has been treated to remove any residual contents of the hazardous substance;
- or for substances that have a health or environmental hazard, or corrosive to metal, the contents of the residue in the package are below the threshold for the substance to be classified as hazardous in the Hazardous Substances (Hazard Classification) Notice 2020.

**Section 14: Transport information**

**IATA**

<b>UN number or ID number</b>	UN1133
<b>UN proper shipping name</b>	Adhesives solution
<b>Transport hazard class(es)</b>	3
<b>Packing group</b>	II
<b>Description</b>	UN1133, Adhesives solution, 3, II
<b>Special Provisions</b>	A3
<b>ERG Code</b>	3L

**IMDG**

<b>UN number or ID number</b>	UN1133
<b>UN proper shipping name</b>	ADHESIVES SOLUTION
<b>Transport hazard class(es)</b>	3
<b>Packing group</b>	II
<b>Description</b>	UN1133, ADHESIVES SOLUTION (n-Heptane), 3, II, (-4°C C.C.), Marine pollutant
<b>EmS-No.</b>	F-E, S-D

**Section 15: Regulatory information**

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

**National regulations**

**EPA New Zealand HSNO approval code or group standard**      HSR002596 - Laboratory Chemicals and Reagent Kits

**National regulations**

There are no applicable tolerable exposure limits or environmental exposure limits according to the EPA Controls for Hazardous Substances

**Certified handlers, tracking and controlled substance license requirements**

Certified handlers are required for some substances. This includes substances requiring a controlled substance license, and most explosives, vertebrates toxic agents, and certain fumigants. Acutely toxic substances which are a Category 1 or 2, such as pesticides also require Certified handlers. Please check the Health and Safety at Work Act 2015 for further information  
Tracking is required for some highly hazardous substances. These substances need to be under the control of an appropriately trained person or appropriately secured. Please check

the Health and Safety at Work Act 2015 for further information  
Controlled substance licenses are required to possess certain explosives, vertebrate toxic agents and fumigants. See Part 7 of the Health and Safety at Work Regulation 2017 for more information

### International Regulations

**The Montreal Protocol on Substances that Deplete the Ozone Layer** Not applicable

**The Stockholm Convention on Persistent Organic Pollutants** Not applicable

**The Rotterdam Convention** Not applicable

### International Inventories

Contact supplier for inventory compliance status

**TSCA**

## Section 16: Other information

<b>Issuing Date</b>	15-Oct-2021
<b>Revision date</b>	04-Dec-2024
<b>Revision Note</b>	Updated format.

### Key or legend to abbreviations and acronyms used in the safety data sheet

#### Legend

ACGIH	American Conference of Governmental Industrial Hygienists
ADN	Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Europe)
ADR	Agreement concerning the International Carriage of Dangerous Goods by Road (Europe)
AIIC	Australian Inventory of Industrial Chemicals
ATE	Acute Toxicity Estimate
ASTM	American Society for the Testing of Materials
bar	Biological Reference Values for Chemical Compounds in the Work Area
BAT	Biological tolerance values for occupational exposure
BEL	Biological exposure limits
bw	Body weight
Ceiling	Maximum limit value
CMR	Carcinogen, Mutagen or Reproductive Toxicant
DOT	Department of Transportation (United States)
DSL	Domestic Substances List (Canada)
EmS	Emergency Schedule
ENCS	Existing and New Chemical Substances (Japan)
EPA	Environmental Protection Agency
GHS	Globally Harmonized System
IARC	International Agency for Research on Cancer
IATA	International Air Transport Association
IBC	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
ICAO	International Civil Aviation Organisation
IECSC	Inventory of Existing Chemical Substances in China
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
ISO	International Organisation for Standardisation

KECI	Korean Existing Chemicals Inventory
LC50	Lethal Concentration to 50% of a test population
LD50	Lethal Dose to 50% of a test population (Median Lethal Dose)
MARPOL	International Convention for the Prevention of Pollution from Ships
n.o.s.	Not Otherwise Specified
NOAEC	No Observed Adverse Effect Concentration
NOAEL	No Observed Adverse Effect Level
NOELR	No Observable Effect Loading Rate
NZIoC	New Zealand Inventory of Chemicals
OECD	Organisation for Economic Cooperation and Development
OEL	Occupational exposure limits
PBT	Persistent, Bioaccumulative and Toxic substance
PICCS	Philippines Inventory of Chemicals and Chemical Substances
PMT	Persistent, Mobile and Toxic
PPE	Personal protective equipment
QSAR	Quantitative Structure Activity Relationship
RID	Agreement concerning the International Carriage of Dangerous Goods by Rail (Europe)
SADT	Self-Accelerating Decomposition Temperature
SAR	Structure-activity relationship
SDS	Safety Data Sheet
SL	Surface Limit
STEL	Short Term Exposure Limit
STOT RE	Specific target organ toxicity - Repeated exposure
STOT SE	Specific target organ toxicity - Single exposure
TCSI	Taiwan Chemical Substance Inventory
TDG	Transport of Dangerous Goods (Canada)
TSCA	Toxic Substances Control Act (United States)
TWA	Time-Weighted Average
UN	United Nations
VOC	Volatile organic compounds
vPvB	Very Persistent and Very Bioaccumulative
vPvM	Very Persistent and Very Mobile
Sen+	Sensitiser
Sk*	Skin designation
**	Hazard Designation

#### Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)  
U.S. Environmental Protection Agency ChemView Database  
European Food Safety Authority (EFSA)  
Environmental Protection Agency  
Acute Exposure Guideline Level(s) (AEGL(s))  
U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act  
U.S. Environmental Protection Agency High Production Volume Chemicals  
Food Research Journal  
Hazardous Substance Database  
International Uniform Chemical Information Database (IUCLID)  
Japan GHS Classification  
Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)  
NIOSH (National Institute for Occupational Safety and Health)  
National Library of Medicine's ChemID Plus (NLM CIP)  
National Library of Medicine's PubMed database (NLM PUBMED)  
U.S. National Toxicology Program (NTP)  
New Zealand's Chemical Classification and Information Database (CCID)  
Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications  
Organisation for Economic Co-operation and Development High Production Volume Chemicals Program  
Organisation for Economic Co-operation and Development Screening Information Data Set  
World Health Organization

**Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of Safety Data Sheet**